



ken-a-vision[®]
KNOWLEDGE THROUGH VISION

**Megapixel Digital
Video Flex[®]
7200UM**

Ken-A-Vision[®] created the USB, High resolution Megapixel Digital Video Flex with a custom CMOS chip, glass lens and Quick Turn[™] Focus. Complete with Applied Vision[™] Software, this powerful analysis tool enhances learning by capturing details you might otherwise miss.



**Additional
Features and Benefits**

- **Digital Signal** For crisp clear images
- **Software** Applied Vision[™] Software
- **Poly-Carbonate Build** Allows for rugged conditions and environment
- **C-mount Lens** Interchangeable lens system accommodates all c-mount lenses.
- **Quick Focus[™]** 1/4 turn moves from 6mm to infinity
- **Kensington Slot** Provides anti-theft control
- **Discovery Scope Kit** To view specimens without a microscope
- **High Resolution** 1280 x 1024 SXGA

7200UM SPECIFICATIONS

Weight	5.95 lbs / 2.70 Kg	Real-Time Video	up to 30 Frames per second
Base Material	Polycarbonate	Magnification	50:1 and Greater
Base Dimension	7"x7"x7" / 18x18x18 cm	Eyepiece Adapters	28mm & Built-in 34.5 mm
Height	30" / 76.2 cm	Connection	Pre-attached 6' (1.82 m) USB USB 2.0
Head Type	Ball and Socket	Software	Applied Vision [™] Software
Lens	8 mm Glass	Options Included	Discovery Scope Kit & Case
Focus	1/4 turn Quick Focus [™]	Warranty	5 Years
Focal distance	6mm to infinity		
Resolution	1280 x 1024 SXGA		

Specifications subject to change without notice



ken-a-vision[®]
KNOWLEDGE THROUGH VISION

(816) 353-4787 Fax (816) 358-5072

email: info@ken-a-vision.com