

# FlexCam

## by Ken-A-Vision

### StudentCam

910-171-120

The versatile, feature-rich StudentCam® lets teachers share detailed documents, images and video and includes advanced lighting controls, such as fluorescent and brightness adjustments. The StudentCam includes a microscope eyepiece adapters.

A built-in microphone captures audio for presentations or videoconferencing. The StudentCam now offers higher resolution, with 410K pixels and more than 480 TV lines. Light sensitivity has also been increased to only 1.5 Lux, and the signal-to-noise ratio is greater than 48 dB.



### Additional Features and Benefits

- Built-in Microphone
- S-Video Connection
- S-Video to RCA Adapter
- Microscope Adapter Kit
- Flickerless Mode
- Adjustable Brightness Control

#### 910-171-120 SPECIFICATIONS

Weight	4.2 lbs / 1.9 kg	Image Type	CCD
Dimension	7" x 7.5" x 22"	Resolution	480 Lines
	18 cm x 19 cm x 56 cm	Real-Time Video	up to 30 Frames per second
Gooseneck Length	20" (50.8 cm)	Output Signal	Composite
Lens	8 mm C-mount		S-Video
Focal Range	1/4" - 50' / .6 cm - 15.24 cm		Audio
Illumination	3 Lux	Microscope Adapters	26 mm, 28 mm & 34 mm
Iris	Electronic	Warranty	2 Years
Grain Control	Auto		

Specifications subject to change without notice



(816) 353-4787 Fax (816) 358-5072

email: [info@ken-a-vision.com](mailto:info@ken-a-vision.com)

kav.fly.StudentCam.905v1