



ken-a-vision[®]
KNOWLEDGE THROUGH VISION

Cordless Dual Purpose Comprehensive Scope T-1931C

The new Cordless Comprehensive Dual Purpose Scope combines a standard compound microscope with limited stereo microscope capability for one low price.



Features and Benefits

- Separate Top Light for Viewing Large, Macro or Opaque Specimens
- Bottom Light for Standard Microscope Slides
- Light Rheostat for Better Light Control
- Monocular Head



T-1931C SPECIFICATIONS

Weight	15 lbs / 8 Kg	Stage	x/y Floating
Base Material	All metal construction	Above Stage Lighting	Yes
Dimension	8.5"L X 6.5"W X 15"H	Below Stage Lighting	Yes
	21.6 cm X 16.5 cm X 38.1 cm	Light Source	Ken-A-Vision Cool Lighting
Head Type	Monocular	Bulb Life	100,000 hours
Eyepiece	10x Widefield	Power Source	Rechargeable NiMH Batteries
Objective Lens	2x, 4x, 10x, 40xR DIN	Warranty	10 Years, Excluding Batteries & Lamp
Focal Adjustment	Coaxial Coarse & Fine		

Specifications subject to change without notice



ken-a-vision[®]
KNOWLEDGE THROUGH VISION

(816) 353-4787 Fax (816) 358-5072

email: info@ken-a-vision.com