



ken-a-vision[®]
KNOWLEDGE THROUGH VISION

Polarizing Comprehensive Scope T-1941

The Comprehensive Polarizing Scope offers high quality construction and a 10 year warranty. This microscope is the ideal microscope for teaching geology, mineralogy, metallurgy, petrology, chemistry or viewing anything enhanced by polarization.



Additional Features and Benefits

- Perfect for Viewing Anisotropic Materials
- Slot for Compensators/Retardation Plates
- 360° Round Rotating Graduation Marked Stage
- Light Rheostat for Better Light Control

T-1941 SERIES SPECIFICATIONS

Weight	11 lbs/ 5 kgs	Analyzer Filter	Yes
Dimension	8.5"L X 6.5"W X 15"H	Coaxial Focal Adjustment	Yes
	21.6 cm X 16.5 cm X 38.1 cm	Iris Diaphragm	Yes
Eyepiece	10x Widefield	Adjustable Abbe Condenser	Yes
Objective Lens	4x, 10x, 40xR DIN	Warranty	10 Years
Head Type	Monocular		
Analyzer Lens in Eyetube	Yes		
Rotating Stage with Clips	Yes		
Polarizer Filter	Yes		

Specifications subject to change without notice



ken-a-vision[®]
KNOWLEDGE THROUGH VISION

(816) 353-4787 Fax (816) 358-5072

email: info@ken-a-vision.com