



**ken-a-vision**<sup>®</sup>  
KNOWLEDGE THROUGH VISION

## Megapixel PupilCAM<sup>®</sup> 1401KEM

Make your microscope digital with the ultra high resolution PupilCAM! The Ken-A-Vision Megapixel PupilCAM 1401KEM is a 1.3 megapixel high resolution microscope camera with digital CMOS technology. Complete with Applied Vision™ Software. It has a factory installed five-foot USB cable and a built-in 10X eyepiece, which you place into the eyepiece tube.



### Features and Benefits

- One Piece Construction
- Plug-N-Play
- Resolution
- Applied Vision™ Software

No parts to lose  
Easy to use  
1280 x 1024 SXGA  
Record images or movies  
Time-lapse sequences or movies  
Still images JPG, BMP, TIFF can be emailed or printed

#### 1401KEM SPECIFICATIONS

Weight	1 lb / 0.45 kgs	Exposed Cable Length	3' 7" (1.12m) USB Cable
Weight Packaged	2 lbs / 0.91 kgs	Resolution	1280 x 1024 SXGA
Image Type	1/3" CMOS	Frame Rate	up to 30 frames per second
Lens	8mm	Software	Applied Vision™
Eyepiece	10X Widefield	Warranty	1 Year
Power Source	From Computer USB Port		
Connection	Factory Installed 5' (1.52 m) USB Cable USB 2.0		

Specifications subject to change without notice



**ken-a-vision**<sup>®</sup>  
KNOWLEDGE THROUGH VISION

(816) 353-4787 Fax (816) 358-5072

email: [info@ken-a-vision.com](mailto:info@ken-a-vision.com)