

# Centrifuges



| Model             | Ultra 8 Variable           | Ultra 8 Select             | Ultra 8 Digital            | E8 Fixed                  | E8 Variable                | Portafuge                 |
|-------------------|----------------------------|----------------------------|----------------------------|---------------------------|----------------------------|---------------------------|
| Model #           | ULC-08AV-1501              | ULC-08AS-1501              | ULC-08AD-1501              | E8C-08AF-1501             | E8C-08AV-1501              | E8C-08AF-150P             |
| Type              | Test Tube                  | Test Tube                  | Test Tube                  | Test Tube                 | Test Tube                  | Test Tube                 |
| Features          | Variable Speed             | 4 Speed Selectable         | Programmable               | Fixed Speed               | Variable Speed             | Fixed Speed               |
| Capacity          | 8 place X 15ml             | 8 place X 15ml             | 8 place X 15ml             | 8 place X 15ml            | 8 place X 15ml             | 8 place X 15ml            |
| Tube Size Min/Max | 10.25 X 65mm / 17 X 130mm  | 10.25 X 65mm / 17 X 130mm  | 10.25 X 65mm / 17 X 130mm  | 10.25 X 65mm / 17 X 130mm | 10.25 X 65mm / 17 X 130mm  | 10.25 X 65mm / 17 X 130mm |
| Max g             | 1400                       | 1400                       | 1400                       | 1400                      | 1400                       | 1500                      |
| Max Speed (rpm)   | 3300                       | 3300                       | 3300                       | 3300                      | 3300                       | +/-3300                   |
| Uses              | Blood, Urine, Semen, Fecal | Blood, Urine, Semen, Fecal | Blood, Urine, Semen, Fecal | Blood Only                | Blood, Urine, Semen, Fecal | Blood Only                |



| Model             | C3 Select                         | C5   | C5 - 24 Place  | CxR                       |
|-------------------|-----------------------------------|--|--|---------------------------|
| Model #           | C3C-08SS-1501                     | C5C-08SV-15T1  | C5C-BKSV-15T3  | CXC-06SF-1501             |
| Type              | Test Tube Swing Out               | Test Tube Swing Out  | Test Tube  | Test Tube Swing Out       |
| Features          | Selectable Pre-set Speeds to 3300 | 5000 RPM Programmable  | 5000 RPM Programmable  | 3400 RPM Short Profile    |
| Capacity          | 8 place X 15ml<br>*4 place X 50ml | 8 place X 15ml<br>*4 place X 50ml  | 24 X 2.5-10ml<br>*4 place X 50ml tube inserts                      | 6 place X 15ml            |
| Tube Size Min/Max | 10.25 X 65mm / 17 X 130mm         | 10.25 X 65mm / 16 X 130mm  | 10.25 X 65mm / 17 X 130mm  | 10.25 X 65mm / 16 X 130mm |
| Max g             | 1500                              | 4000   | 4000   | 1800                      |
| Max Speed (rpm)   | 3000                              | 5000   | 5000   | 3400                      |
| Uses              | Blood, Urine, Semen, Fecal        | Blood (STAT), Urine, Semen, Fecal, Hematocrits, PRP/PPP, Stem-Cell, Irap | Blood (STAT), Urine, Semen, Fecal, Hematocrits, PRP/PPP, Stem-Cell | Blood Only                |



\*Optional Features

P 770.270.1394  
F 770.270.2389

865 Marathon Parkway  
Lawrenceville GA 30046

LWScientific.com

# Centrifuges



| Model             | V24 Combo   | V24 Combo Tach  | M24 Hematocrit        | Zipocrit                   | ZipSpin                      |
|-------------------|---|---|-----------------------|----------------------------|------------------------------|
| Model #           | CMC-24AV-1501<br>(choose rotor)                         | CMC-24AV-15T1<br>(choose rotor)                         | HEC-24HF-7501         | ZPC-04HF-750P              | ZPC-06AF-020P                |
| Type              | Combination   | Combination   | Microhematocrit       | Microhematocrit            | Microcentrifuge              |
| Features          | Test Tube & Microhematocrit                             | Test Tube & Microhematocrit                             | Heavy Duty 24 Place   | Portable angled hematocrit | Portable 6 & 12 place rotors |
| Capacity          | 24 X 15ml or 24 X 75mm                                  | 24 X 15ml or 24 X 75mm                                  | 24 X 75mm             | 4 X 75mm                   | (6 & 12) X 2ml               |
| Tube Size Min/Max | 10.25 X 65mm /17X130mm or 75mm capillary                | 10.25X65mm /17X130mm or 75mm capillary                  | 75mm capillary (only) | 75mm capillary (only)      | 0.5ml / 2.0ml                |
| Max g             | 2450 - 14900  | 2450 - 14900  | 14900                 | 4400                       | 2800                         |
| Max Speed (rpm)   | 4000 - 12000  | 4000 - 12000  | 4000 - 12000          | 12000                      | 11000                        |
| Users             | Blood, Urine, Semen, Fecal, Hematocrits, DNA Separation | Blood, Urine, Semen, Fecal, Hematocrits, DNA Separation | Hematocrits Only      | Hematocrits Only           | Blood, Urine, DNA Separation |



P 770.270.1394  
F 770.270.2389

865 Marathon Parkway  
Lawrenceville GA 30046

LWScientific.com