



**ken-a-vision®**  
KNOWLEDGE THROUGH VISION

## Comprehensive Scope 2 T-19011C

The Comprehensive Scope 2 series incorporates a professional instrument priced within an educational budget. This unit is suitable for ages 12 and up including university level courses and has the advantage of Ken-A-Vision's own superior cordless technology. These units have a modern, modular design on a very stable base. The arm of the microscope includes a built in handle for ease of carrying and an integrated storage system for the power cord. Accessory kits make available 14 different reconfigurations of microscope features.



### Features and Benefits

- ◆ 10x Wide Field Eyepiece
- ◆ Achromat Objectives
- ◆ Dimmer (Rheostat) Switch
- ◆ Reverse Nosepiece
- ◆ Coaxial Focusing Knob
- ◆ Carrying Handle
- ◆ Cordless
- ◆ Rechargeable Batteries
- ◆ Storage Door
- ◆ Ken-A-Vision Cool Lighting
- 19mm wide eyepiece
- 4x, 10x, 40xR DIN
- Easy adjustment of transmitted light
- Ease of slide insertion and manipulation
- Separate coarse and fine adjusting
- Convenience
- Use with or without charger
- NiMH batteries recharge in 8, last up to 40 hours
- Charging adapter stays with unit
- Patented LED bright white light, long bulb life

### T-19011C SPECIFICATIONS

Weight	10.4lbs/ 4.7kgs	Objective Lens	4x, 10x, 40xR DIN
Size	10" x 7.13" x 14"	Stage	Stage Clips
	25.4cm x 18cm x 35.6cm	Focus	Separate Coarse & Fine
Head Type	Monocular	Diaphragm	5 Holed Disk
Eyepiece	10x Widefield with Pointer	Light Source	Ken-A-Vision Cool Lighting
Nosepiece	Reversed Four Hole	Variable Light Control	Yes - Rheostat
		Power	Cordless - NiMH Batteries
		Warranty	10 years

Specifications subject to change without notice



**ken-a-vision®**  
KNOWLEDGE THROUGH VISION

(816) 353-4787 Fax (816) 358-5072

email: [info@ken-a-vision.com](mailto:info@ken-a-vision.com)